



Solutions Catalog

HOSPITALS & SURGICAL CENTERS

www.ansell.com

Ansell



Ansell is an innovator in high-quality, high-performance, custom healthcare safety solutions. For more than a century, we’ve been setting industry milestones with glove and safety device advances that consistently provide for you and your patients’ safety, well-being and peace of mind.

TABLE OF CONTENTS

Our Promise	3
Our History	5
Our Solutions	7
Surgical Gloves	9
Examination Gloves	19
Safety Solutions	25
AnsellCARES	57



OUR PROMISE

PRODUCT INNOVATION

Solutions for today and tomorrow

ADVANCED TECHNOLOGIES

Ongoing research & development delivers ground-breaking advances in materials, product design and manufacturing processes.

BEST PRACTICES

Continuous investment in manufacturing, quality assurance and supply chain so we can deliver superior products when you need them.

SUPPLY MANAGEMENT

Helping you contain costs by simplifying product selection and streamlining inventory while continuing to meet all caregivers' needs.

MARKET FOCUS

Connected to caregivers and patients

CONSULTANT NETWORK

Global network of Key Opinion Leaders helps ensure reliable products that reduce infection risks & promote safety and quality outcomes.

PRODUCT & CLINICAL EXPERTISE

Experienced sales & clinical specialists stay abreast of healthcare technologies and best practices to provide you with evidence-based safety recommendations.

CONTINUING EDUCATION

Accredited CE programs on caregiver and patient safety topics and other free learning resources available anytime on our website.



GLOBAL LEADERSHIP

Comprehensive customer support

WORLDWIDE DISTRIBUTION

Reliable network of international and local distributors deliver Ansell protective solutions directly to you.

BRAND LOYALTY

A global brand you can trust — with #1 global market share for surgical gloves.

GIVING BACK

From product donations and volunteerism to energy efficiency, we're committed to engaging with the communities we serve and leading our industry by example.

OUR SOLUTIONS

In our pursuit of developing solutions for healthcare safety, we’ve always looked beyond what is currently feasible and set our sights on what should be feasible. The results are superior gloves and safety devices with ground breaking, proprietary technologies, some of which have been patented, that offer you and your patients optimal protection in any healthcare environment and for any user need.

ANSELL PROPRIETARY TECHNOLOGIES

 **A.R.T.**
Pinhole Barrier Technology

This technology is proprietary process, integrated into the manufacturing line that helps add an extra level of robustness by minimizing possible defects that may later result in pinholes.

 **AMT**
Antimicrobial Technology

This technology enables the use of chlorhexidine gluconate (CHG) on the inside of the glove to suppress the regrowth of bacteria on hands and prevent bacterial and viral migration in the event of a glove breach.

 **DERMASHIELD**
Donning Aid Technology

This technology is a proprietary inner coating for fast and easy damp or dry hand donning whether it is the first don or during intraoperative changes.

 **ERGOFORM**
Ergonomic Design Technology

This technology enables Ansell to design safety solutions that support musculoskeletal health during repetitive tasks to improve worker performance.

 **HYDRASOFT**
Skin Moisturizing Technology

This technology is an odorless, water-based coating applied on the inside of the glove that helps retain moisture and rehydrate dry skin to improve the wearers’ skin condition.

 **PI-KARE**
Skin-friendly PI Technology

This technology enables the elimination of standard chemical accelerators from polyisoprene (PI) gloves known to cause chemical Type IV allergies and sensitivities, making them safer to use.

 **SENSOPRENE**
Ultra Thin Neoprene Technology

This technology delivers a glove up to 30% thinner than Ansell’s standard neoprene surgical gloves so wearers’ experience a high degree of sensitivity without compromising on durability and the optimum skin protection needed from allergies and sensitivities caused by latex or chemical accelerators.

SURGICAL GLOVES

GAMMEX® ENCORE®

Driven by innovation, Ansell’s range of surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



EXAMINATION GLOVES

MICROTOUCH®

Ansell’s premium examination gloves are scientifically developed to fully equip and comfortably protect medical professionals so that they have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



SAFETY SOLUTIONS

SANDEL®

Developed in collaboration with clinical professionals, SANDEL products are designed to address key high-risk areas in the healthcare environment. Our signature bright orange products help reduce the risk of preventable injuries, errors and infections. SANDEL products equip caregivers with solutions to do their job safely while improving patient outcomes.



GAMMEX®
ENCORE®

SURGICAL GLOVES

LOW ALLERGENICITY

Skin-friendly polyisoprene (PI) gloves free from chemical accelerators known to cause chemical allergies, 100% chemical accelerator-free neoprene synthetics and low-protein latex range minimize your risk of developing sensitivities and triggering allergic reactions that can hinder your performance.

HAND HEALTH

Exclusive formers ergonomically engineered to custom-fit the natural shape of your hands reduce the risk of repetitive stress and muscle fatigue. Innovative, proprietary formulations are softer and thinner yet remarkably strong, optimizing comfort without compromising barrier effectiveness.

USER-FOCUSED TECHNOLOGIES

Proprietary technologies are designed to enrich your user experience and enhance performance—including skin moisturizing, damp donning, secure cuffs, enhanced grip, double gloving systems and more.

ECO-FRIENDLY PACKAGING

Patented SMART Pack meets the priorities of today's busy O.R. environment: space-saving and versatile design adaptable to any environment, smaller size and lightweight, prominent labeling and color coding for easy identification, effortless opening with easy return of unused gloves, and less plastic and cardboard waste.

ALLERGY & MATERIAL GUIDE

Please note that the following descriptions and tables are representative of basic material characteristics. Individual glove performance may vary.

	PI	NEOPRENE	LATEX
LEVEL OF BARRIER PROTECTION	Excellent	Excellent	Excellent
ALLERGEN CONTENT	Very Good	Very Good to Excellent	Varies
STRENGTH & DURABILITY	Very Good	Excellent	Excellent
ELASTICITY	Excellent	Excellent	Excellent
FIT, COMFORT	Excellent	Very Good	Excellent
ECONOMY	Fair	Good	Very Good

SYNTHETIC RANGE



	GAMMEX® Non-Latex PI			GAMMEX® Non-Latex PI Green			GAMMEX® Non-Latex PI Textured			GAMMEX® Non-Latex PI Underglove			GAMMEX® Non-Latex PI Micro			GAMMEX® Non-Latex PI Ortho		
MATERIAL	PI (Polyisoprene)			PI (Polyisoprene)			PI (Polyisoprene)			PI (Polyisoprene)			PI (Polyisoprene)			PI (Polyisoprene)		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none"> • Softer and ‘stretchier’ PI glove featuring superior comfort compared to Ansell’s latex gloves • Exceptional instrument handling and grip • Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities for optimized allergenic profile • Effortless dry and damp donning as well as intra-operative changes 			<ul style="list-style-type: none"> • Softer and ‘stretchier’ PI glove featuring superior comfort compared to Ansell’s latex gloves • Optimal instrument handling and grip • Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities for optimized allergenic profile • Tested for use with chemotherapy drugs and a wide range of chemicals 			<ul style="list-style-type: none"> • Exceptional grip under dry and damp conditions • Softer and ‘stretchier’ PI glove featuring superior comfort compared to Ansell’s latex gloves • Ideal to be used as an outer glove when double gloving • Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities • Effortless dry and damp donning as well as intra-operative changes 			<ul style="list-style-type: none"> • Optimal comfort when double gloving • Green color for easy glove breach identification • Softer and ‘stretchier’ PI glove featuring superior comfort compared to Ansell’s latex gloves • Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities • It can be used as a single glove 			<ul style="list-style-type: none"> • 15% thinner on average than Ansell’s standard polyisoprene surgical gloves for unsurpassed dexterity and tactile sensitivity • Softer and ‘stretchier’ PI glove featuring superior comfort compared to Ansell’s latex gloves • Exceptional instrument handling and grip • Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities 			<ul style="list-style-type: none"> • 50% thicker on average than Ansell’s standard powder-free latex surgical gloves • Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities • Unique surface treatment for exceptional handling and grip • Effortless dry and damp donning as well as intra-operative changes 		
STERILE/NON-STERILE	Sterile			Sterile			Sterile			Sterile			Sterile			Sterile		
POWDER CONTENT	Powder-free			Powder-free			Powder-free			Powder-free			Powder-free			Powder-free		
GLOVE COLOR	White			Green			White			Green			White			Green		
INTERNAL GLOVE SURFACE	Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology		
EXTERNAL GLOVE SURFACE	Smooth with micro-textured finish			Smooth with micro-textured finish			Textured fingers and palm			Smooth with micro-textured finish			Smooth with micro-textured finish			Smooth with micro-textured finish		
THICKNESS	Standard			Standard			Standard			Micro			Micro			Ortho		
GRIP LEVEL	Moderate			Moderate			High			Low			Moderate			Moderate		
VULCANIZATION CHEMICAL ACCELERATORS	Xanthates, Zinc dialkyldithiocarbamate			Xanthates, Zinc dialkyldithiocarbamate			Xanthates, Zinc dialkyldithiocarbamate			Xanthates, Zinc dialkyldithiocarbamate			Xanthates, Zinc dialkyldithiocarbamate			Xanthates, Zinc dialkyldithiocarbamate		
ALLERGY PREVENTION	Latex (Type I)			Latex (Type I)			Latex (Type I)			Latex (Type I)			Latex (Type I)			Latex (Type I)		
DOUBLE GLOVING RECOMMENDATION	Outer glove or underglove			Outer glove or underglove			Outer glove			Underglove			Outer glove			No		
FREEDOM FROM HOLES (AQL)	Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)		
ANSELL TECHNOLOGIES	PI-KARE™, ERGOFORM™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			PI-KARE™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			PI-KARE™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			PI-KARE™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			PI-KARE™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			PI-KARE™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	5.5	20685755	50 pairs per box;	5.5	20685255	50 pairs per box;	5.5	20688255	50 pairs per box;	5.5	20687255	50 pairs per box;	5.5	20685955	50 pairs per box;	5.5	N/A	50 pairs per box;
	6	20685760		6	20685260		6	20688260		6	20687260		6	20685960		6	20686560	
	6.5	20685765		6.5	20685265		6.5	20688265		6.5	20687265		6.5	20685965		6.5	20686565	
	7	20685770	4 boxes per carton/case;	7	20685270	4 boxes per carton/case;	7	20688270	4 boxes per carton/case;	7	20687270	4 boxes per carton/case;	7	20685970	4 boxes per carton/case;	7	20686570	4 boxes per carton/case;
	7.5	20685775		7.5	20685275		7.5	20688275		7.5	20687275		7.5	20685975		7.5	20686575	
	8	20685780	200 pairs per carton/case	8	20685280	200 pairs per carton/case	8	20688280	200 pairs per carton/case	8	20687280	200 pairs per carton/case	8	20685980	200 pairs per carton/case	8	20686580	200 pairs per carton/case
	8.5	20685785		8.5	20685285		8.5	20688285		8.5	20687285		8.5	20685985		8.5	20686585	
	9	20685790		9	20685290		9	20688290		9	20687290		9	20685990		9	20686590	

* Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

SYNTHETIC RANGE



	GAMMEX® PI Hybrid			GAMMEX® PI Plus Glove-in-Glove™ System			GAMMEX® Non-Latex			GAMMEX® Non-Latex Sensitive		
MATERIAL	PI (Polyisoprene), Neoprene (Polychloroprene)			Outer: PI (Polyisoprene) ; GAMMEX® Non-Latex PI Textured Inner: PI (Polyisoprene) ; GAMMEX® Non-Latex PI Underglove			Neoprene (Polychloroprene)			Neoprene (Polychloroprene)		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none"> Patent-pending HYBRID™ Technology delivers latex-like comfort and high durability Ideal for use as an outer glove when double gloving Non-latex to eliminate the risk of latex Type I allergy Zinc 2-mercaptobenzothiazole-free (ZMBT-free), diphenylguanidinefree (DPG-free) and cetylpyridinium chloride-free (CPC-free) 			<ul style="list-style-type: none"> Exceptional grip under dry and damp conditions Softer and 'stretchier' PI glove featuring superior comfort compared to Ansell's latex gloves Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities Outer and inner gloves are pre-donned and aligned at the finger tips for easy double gloving in one don, saving on half the donning time 			<ul style="list-style-type: none"> Non-latex and 100% chemical accelerator-free Extended cuff length for extra protection Optimal instrument handling and grip Tested for use with chemotherapy drugs and a wide range of chemicals* 			<ul style="list-style-type: none"> Non-latex and 100% chemical accelerator-free Featuring SENSOPRENE™ Formulation to deliver a thinner yet robust glove for added comfort and durability Up to 30% thinner on average than Ansell's standard neoprene surgical gloves Tested for use with chemotherapy drugs and a wide range of chemicals* 		
STERILE/NON-STERILE	Sterile			Sterile			Sterile			Sterile		
POWDER CONTENT	Powder-free			Powder-free			Powder-free			Powder-free		
GLOVE COLOR	Natural			Outer: Natural, Inner: Green			Green			White		
INTERNAL GLOVE SURFACE	Polymer Coated with DERMASHIELD™ Technology			Polymer Coated with DERMASHIELD™ Technology			Polymer Coated with DERMASHIELD™ Technology			Polymer Coated with DERMASHIELD™ Technology		
EXTERNAL GLOVE SURFACE	Smooth with micro-textured finish			Smooth with micro-textured finish			Smooth with micro-textured finish			Smooth with micro-textured finish		
THICKNESS	Standard			Outer: Standard; Inner: Micro			Standard			Micro		
GRIP LEVEL	Moderate			Outer: Moderate; Inner: Low			Moderate			Moderate		
VULCANIZATION CHEMICAL ACCELERATORS	Zinc diethyldithiocarbamate (ZDEC), Zinc dibutyldithiocarbamate (ZDBC), Sodium di-n-butylthiocarbamate (SDBC), Xanthogen polysulphides			Xanthates, Zinc dialkyldithiocarbamate			None			None		
ALLERGY PREVENTION	Latex (Type 1)			Latex (Type 1)			Latex (Type 1) & Chemical (Type IV)			Latex (Type 1) & Chemical (Type IV)		
DOUBLE GLOVING RECOMMENDATION	Outer glove			N/A			Outer glove or underglove			Outer glove		
FREEDOM FROM HOLES (AQL)	Manufacturing final release: 0.65 (Inspection level GII)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared		
ANSELL TECHNOLOGIES	HYBRID™, ERGOFORM™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			PI-KARE™, ERGOFORM™, DERMASHIELD™, A.R.T™, GLOVE-IN-GLOVE™, PEARL™, SUREFIT™			DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™			SENSOPRENE™, ERGOFORM™, DERMASHIELD™, A.R.T™, PEARL™, SUREFIT™		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	5.5	340063055	50 pairs per box;	5.5	340082055	2 pairs per pouch;	5.5	8511	50 pairs per box;	5.5	20277255	50 pairs per box;
	6	340063060		6	340082060		6	8512		6	20277260	
	6.5	340063065		6.5	340082065		6.5	8513		6.5	20277265	
	7	340063070	4 boxes per carton/case;	7	340082070	4 boxes per carton/case;	7	8514	4 boxes per carton/case;	7	20277270	4 boxes per carton/case;
	7.5	340063075		7.5	340082075		7.5	8515		7.5	20277275	
	8	340063080	200 pairs per carton/case	8	340082080	200 pairs per carton/case (100 pouches)	8	8516	200 pairs per carton/case	8	20277280	200 pairs per carton/case
	8.5	340063085		8.5	340082085		8.5	8517		8.5	20277285	
	9	340063090		9	340082090		9	8518		9	20277290	

* Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

LATEX RANGE



	ENCORE® Latex Acclaim®			ENCORE® Latex Micro			ENCORE® Latex Moisturizing			ENCORE® Latex Ortho			ENCORE® Latex Texture®			ENCORE® Latex Underglove		
MATERIAL	Natural Rubber Latex			Natural Rubber Latex			Natural Rubber Latex			Natural Rubber Latex			Natural Rubber Latex			Natural Rubber Latex		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none">• Soft and sensitive yet robust surgical glove• Unique surface treatment for exceptional instrument handling and grip• Effortless dry and damp donning as well as intra-operative changes• Tested for use with chemotherapy drugs*			<ul style="list-style-type: none">• 20% thinner on average than Ansell’s standard latex powder-free surgical gloves for unsurpassed dexterity and tactile sensitivity• Its brown color helps reduce glare of surgical lights and makes it suitable to be used as an underglove• The DERMASHIELD™ Technology facilitates dry or damp donning and intra-operative changes			<ul style="list-style-type: none">• Helps keep skin moist and supple• Optimal tactile sensitivity and protection• Soft and sensitive yet a robust surgical glove• Effortless dry and damp donning as well as intra-operative changes			<ul style="list-style-type: none">• 50% thicker on average than Ansell’s standard surgical gloves• Superior grip for dry and damp applications• Effortless dry and damp donning as well as intra-operative changes			<ul style="list-style-type: none">• Exceptional grip under dry and damp conditions• Optimal tactile sensitivity and protection• Effortless dry and damp donning as well as intra-operative changes			<ul style="list-style-type: none">• Helps keep skin moist and supple• Effortless donning and optimal comfort when double gloved• Easy identification of a glove breach• Effortless dry and damp donning as well as intra-operative changes		
STERILE/NON-STERILE	Sterile			Sterile			Sterile			Sterile			Sterile			Sterile		
POWDER CONTENT	Powder-free			Powder-free			Powder-free			Powder-free			Powder-free			Powder-free		
GLOVE COLOR	Natural			Brown			Natural			Brown			Natural			Green		
INTERNAL GLOVE SURFACE	Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology			Coated with HYDRASOFT™ Technology			Polymer coated with DERMASHIELD™ Technology			Polymer coated with DERMASHIELD™ Technology			Coated with HYDRASOFT™ Technology		
EXTERNAL GLOVE SURFACE	Smooth			Smooth with micro-textured finish			Smooth			Textured fingers and palm			Textured fingers and palm			Smooth with micro-textured finish		
THICKNESS	Standard			Micro			Standard			Ortho			Standard			Micro		
GRIP LEVEL	Moderate			Moderate			Moderate			High			High			Low		
VULCANIZATION CHEMICAL ACCELERATORS	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)			Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)			Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)			Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)			Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)			Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)		
ALLERGY PREVENTION	None			None			None			None			None			None		
DOUBLE GLOVING RECOMMENDATION	Outer glove			Outer glove or underglove			Outer glove			Underglove			Outer glove			Underglove		
FREEDOM FROM HOLES (AQL)	Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			No, not tested in accordance with ASTM D6978			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)		
ANSELL TECHNOLOGIES	DERMASHIELD™, A.R.T™, SUREFIT™			DERMASHIELD™, A.R.T™, SUREFIT™			HYDRASOFT™, DERMASHIELD™, A.R.T™, SUREFIT™			DERMASHIELD™, A.R.T™, SUREFIT™			DERMASHIELD™, A.R.T™, SUREFIT™			HYDRASOFT™, DERMASHIELD™, A.R.T™, SUREFIT™		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	5.5	5795000	50 pairs per box;	5.5	5787000	50 pairs per box;	5.5	2018655	50 pairs per box;	5.5	N/A	50 pairs per box;	5.5	5785000	50 pairs per box;	5.5	2018455	50 pairs per box;
	6	5795001		6	5787001		6	2018660		6	5788001		6	5785001		6	2018460	
	6.5	5795002		6.5	5787002		6.5	2018665		6.5	5788002		6.5	5785002		6.5	2018465	
	7	5795003	4 boxes per carton/case;	7	5787003	4 boxes per carton/case;	7	2018670	4 boxes per carton/case;	7	5788003	4 boxes per carton/case;	7	5785003	4 boxes per carton/case;	7	2018470	4 boxes per carton/case;
	7.5	5795004		7.5	5787004		7.5	2018675		7.5	5788004		7.5	5785004		7.5	2018475	
	8	5795005	200 pairs per carton/case	8	5787005	200 pairs per carton/case	8	2018680	200 pairs per carton/case	8	5788005	200 pairs per carton/case	8	5785005	200 pairs per carton/case	8	2018480	200 pairs per carton/case
	8.5	5795006		8.5	5787006		8.5	2018685		8.5	5788006		8.5	5785006		8.5	2018485	
	9	5795007		9	5787007		9	2018690		9	5788007		9	5785007		9	2018490	

* Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

LATEX RANGE



	ENCORE® Perry® Style 42® PF			ENCORE® Sensi-Touch® PF		
MATERIAL	Natural Rubber Latex			Natural Rubber Latex		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none">• Ideal comfort and fit• Polymer coated for fast and easy donning• Optimized balance of protection and tactile sensitivity			<ul style="list-style-type: none">• Ideal comfort and fit• Polymer coated for fast and easy donning• Features a smooth finish with the addition of a beaded cuff to help prevent cuff roll-down		
STERILE/NON-STERILE	Sterile			Sterile		
POWDER CONTENT	Powder-free			Powder-free		
GLOVE COLOR	Natural			Natural		
INTERNAL GLOVE SURFACE	Polymer coated			Polymer coated		
EXTERNAL GLOVE SURFACE	Smooth with micro-textured finish			Smooth with micro-textured finish		
THICKNESS	Standard			Standard		
GRIP LEVEL	Moderate			Moderate		
VULCANIZATION CHEMICAL ACCELERATORS	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)			Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)		
ALLERGY PREVENTION	None			None		
DOUBLE GLOVING RECOMMENDATION	Outer glove			Outer glove		
FREEDOM FROM HOLES (AQL)	Manufacturing final release: 0.65 (Inspection level GI)			Manufacturing final release: 0.65 (Inspection level GI)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	No, not tested in accordance with ASTM D6978			No, not tested in accordance with ASTM D6978		
ANSELL TECHNOLOGIES	None			None		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	5.5	5711100PF	50 pairs per box;	5.5	7821PF	50 pairs per box;
	6	5711101PF		6	7822PF	
	6.5	5711102PF		6.5	7823PF	
	7	5711103PF	4 boxes per carton/case;	7	7824PF	4 boxes per carton/case;
	7.5	5711104PF		7.5	7825PF	
	8	5711105PF		8	7826PF	
	8.5	5711106PF	200 pairs per carton/case	8.5	7827PF	200 pairs per carton/case
	9	5711107PF		9	7828PF	

CAUTION: This product contains natural rubber latex which may cause allergic reactions.

SPECIALTY RANGE



	GAMMEX® PI Radiation Attenuation			ENCORE® Radiation Attenuation*			GAMMEX® Cut-Resistant Glove Liner		
MATERIAL	PI (Polyisoprene) with Bismuth Oxide			Natural Rubber Latex with Bismuth Oxide			High Density Fiber: DYNEEMA® DIAMOND, Nylon, Spandex & Polyester		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none">• Reduces radiation skin-dose exposure by up to 53.4% (beam energy level 60 kVp) versus 1% when double gloving with conventional surgical gloves²• Eliminates the risk of developing Type I latex allergy and sensitization• Optimized balance of protection and comfortable fit• Micro-textured finish for good grip control and handling			<ul style="list-style-type: none">• Reduces radiation skin-dose exposure by up to 53.7% (beam energy level 60 kVp) verses 1% when double gloving with conventional surgical gloves²• Optimised protection, comfortable fit and tactile sensitivity• Micro-textured finish for better grip control and handling			<ul style="list-style-type: none">• Adequate cut-resistance and protection during rigorous procedures**• Sterile• ANSI Cut Level 2 as tested by ASTM F1790• Single use		
STERILE/NON-STERILE	Sterile			Sterile			Sterile		
POWDER CONTENT	Powder-free			Powder-free			Powder-free		
GLOVE COLOR	Black			Khaki			White		
INTERNAL GLOVE SURFACE	Chlorinated			Polymer Coated			None		
EXTERNAL GLOVE SURFACE	Textured fingers			Smooth with micro-textured finish			Knitted fabric		
THICKNESS	Standard			Standard			Micro		
GRIP LEVEL	Moderate			Moderate			Low		
VULCANIZATION CHEMICAL ACCELERATORS	Proprietary Information - Not Available			Proprietary Information - Not Available			Proprietary Information - Not Available		
ALLERGY PREVENTION	Latex (Type I)			None			Latex (Type I) & Chemical (Type IV)		
DOUBLE GLOVING RECOMMENDATION	Underglove			Underglove			Underglove		
FREEDOM FROM HOLES (AQL)	Manufacturing final release: 1.5 (Inspection level GI)			Manufacturing final release: 1.0 (Inspection level GI)			• Manufacturing final release: 1.0 (Inspection level GI)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	No, not tested in accordance with ASTM D6978			No, not tested in accordance with ASTM D6978			Not applicable		
ANSELL TECHNOLOGIES	PI-KARE™			None			None		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	5.5	N/A	5 pairs per box	5.5	N/A	5 pairs per box	XS/5-5.5	5789910	5 pairs per box
	6	N/A		6	N/A		S/6-6.5	5789911	
	6.5	20873565		6.5	20873065		M/7-7.5	5789912	
	7	20873570		7	20873070		L/8-8.5	5789913	
	7.5	20873575		7.5	20873075		XL/9	5789914	
	8	20873580		8	20873080				
	8.5	20873585		8.5	20873085				
	9	20873590		9	20873090				

1. Wagner LK and Mulhern OR. Radiation-attenuating surgical gloves: effects of scatter and secondary electron production. Radiology 1996; 200:45–48.

*CAUTION: Radiation Attenuation gloves are designed to protect hands from scattered radiation exposure originating from the X-ray beam during fluoroscopic procedures. They are not intended for use in direct or primary X-ray beams.

**Sharps Injuries in the Operating Room; Massachusetts Sharps Injury Surveillance Data, 2004. Occupational Health Surveillance Program, Massachusetts Department of Public Health. April, 2008
IMPORTANT: Classified as a personal protective shield and are not a medical device. The gloves must be used in combination with medical gloves if they come in contact with blood or other bodily fluids.

WARNING: No liner provides complete protection against cuts. Use care when handling sharp objects.

ERGONOMIC DESIGN

Uniquely engineered using exclusive formers that create a custom-fit with superior comfort; proven to minimize repetitive motion stress and hand muscle fatigue,* helping you avoid pain and injuries.

GRIP & SENSITIVITY

Innovative formulations that are softer and thinner preserve your tactile sensitivity and dexterity; improved grip in both wet and dry conditions helps you work efficiently with more precision.

HAND HEALTH

Ongoing technological advancements ensure you will find a hand protection solution that meets your needs — low-allergen, powder-free latex; accelerator-free synthetics; skin moisturizing coatings and more.

BARRIER PROTECTION

Breakthrough glove materials are lighter and thinner yet meet the most stringent quality benchmarks so you can be confident that you're protecting your patients and yourself.

ALLERGY & MATERIAL GUIDE

Please note that the following descriptions and tables are representative of basic material characteristics. Individual glove performance may vary.

	NITRILE	NEOPRENE	LATEX	VINYL
LEVEL OF BARRIER PROTECTION	Excellent	Excellent	Excellent	Poor
ALLERGEN CONTENT	Very Good	Very Good to Excellent	Varies	Good
STRENGTH & DURABILITY	Excellent	Very Good	Excellent	Poor
ELASTICITY	Very Good	Excellent	Excellent	Fair to Poor
FIT, COMFORT	Very Good	Very Good	Excellent	Fair
ECONOMY	Very Good	Good	Very Good	Very Good

* Data on file.


SYNTHETIC RANGE



	MICROTOUCH® Nitrile			MICROTOUCH® NitraFree™			MICROTOUCH® NitraFree™ Blue			MICROTOUCH® Micro-Thin™ Nitrile			MICROTOUCH® Nitrile E.P.®			MICROTOUCH® Nitrile Sterile Pairs			MICROTOUCH® Central Supply - SPD Plus			MICROTOUCH® Affinity™		
MATERIAL	Nitrile			Nitrile			Nitrile			Nitrile			Nitrile			Nitrile			Nitrile			Neoprene (Polychloroprene)		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none">• Exceptional comfort and fit• Tested for use with over 55 chemicals and approved for use with chemotherapy drugs*• A better than industry standards AQL delivers consistent, high quality production			<ul style="list-style-type: none">• Proudly supports The Breast Cancer Research Foundation®• Protects from latex Type I allergy• Minimizes chemical Type IV allergy• Tested for use with chemotherapy drugs*			<ul style="list-style-type: none">• Manufactured without sulfur-based chemical accelerators, helping to protect users from Type I latex allergies and Type IV chemical allergies• Thin and light weight formulation provides enhanced tactile sensitivity and comfort for longer wear time• Unique berry blue color provides contrast, while masking blood and other stains• Tested for use with chemotherapy drugs			<ul style="list-style-type: none">• Excellent tactile sensitivity is ideal for palpating veins and starting IVs• The textured surface of the glove promotes sure handling of materials and instruments• 300 gloves per box reduces restocking and reorders			<ul style="list-style-type: none">• Extended 12” cuff provides added protection• Excellent chemical and puncture resistance• Tested for use with chemotherapy drugs*			<ul style="list-style-type: none">• Packed as a pair of gloves per pack• 12” extended cuff for additional protection• Gamma sterilized minimum 25kGy to provide sterility• Tested for use with chemotherapy drugs*			<ul style="list-style-type: none">• 16in / 400mm long to extend protection to the elbow• Tested against 10 commonly used chemicals in the Central Supply/SPD Department• Robust 9mil thickness at the finger for enhanced protection against tears and puncture• Disposable to minimize cross-contamination occurrences• Offer cost savings and improved protection versus fluid resistant cuff sleeves			<ul style="list-style-type: none">• Robust Neoprene formulation provides comfort and excellent tear and puncture resistance• Textured fingers for secure grip• Protects from Type I latex allergy• Tested for use with chemotherapy drugs*		
STERILE/NON-STERILE	Non-Sterile			Non-Sterile			Non-Sterile			Non-Sterile			Non-Sterile			Sterile			Non-Sterile			Non-Sterile		
GLOVE COLOR	Blue			Pink			Blue			Blue			Blue			Blue			Blue			Green		
CUFF LENGTH	Standard			Standard			Standard			Extended			Extended			Extended			Long			Standard		
POWDER CONTENT	Powder-free			Powder-free			Powder-free			Powder-free			Powder-free			Powder-free			Powder-free			Powder-free		
EXTERNAL GLOVE SURFACE	Textured fingers			Textured fingers			Textured fingers			Textured fingers			Textured fingers			Textured fingers			Textured fingers			Textured fingers		
FREEDOM FROM HOLES (AQL)	1.5			1.5			1.5			1.5			1.5			1.5			1.5			1.5		
CUFF THICKNESS	0.05 mm			0.07 mm			0.07 mm			0.05 mm			0.08 mm			0.08 mm			0.10 mm			0.11 mm		
FINGER THICKNESS	0.09 mm			0.14 mm			0.10 mm			0.09 mm			0.16 mm			0.16 mm			0.24 mm			0.16 mm		
PALM THICKNESS	0.07 mm			0.08 mm			0.07 mm			0.06 mm			0.11 mm			0.11 mm			0.21 mm			0.13 mm		
GLOVE LENGTH	240 mm			240 mm			240 mm			230 mm			295 mm			295 mm			400 mm			240 mm		
ALLERGY PREVENTION	Latex (Type 1)			Latex (Type 1) & Chemical (Type IV)			Latex (Type 1) & Chemical (Type IV)			Latex (Type 1)			Latex (Type 1)			Latex (Type 1)			Latex (Type 1)			Latex (Type 1)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 and US FDA cleared			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)			Yes, in accordance with ASTM D6978 and US FDA cleared		
ANSELL TECHNOLOGIES	N/A			N/A			N/A			N/A			N/A			N/A			N/A			N/A		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	XS	6034300	200 gloves per box (size XS-XL)	XS	6034510	100 gloves per box	XS	313028060	200 gloves per box (size XS-L)	XS	6034310	300 gloves per box	XS	6034050	100 gloves per box	XS	N/A	50 pairs per box	S	ULNCS-1	50 gloves per box	XS	3770	100 gloves per box
	S	6034301		S	6034511		S	313028070		S	6034311		S	6034051		S	6034151		M	ULNCS-2		S	3771	
	M	6034302	180 gloves per box (size XXL)	M	6034512	10 boxes per case/carton	M	313028080	180 gloves per box (size XL)	M	6034312	10 boxes per carton/case	M	6034052	10 boxes per carton/case	M	6034152	4 boxes per carton/case	L	ULNCS-3	10 boxes per carton/case	M	3772	10 boxes per carton/case
	L	6034303		L	6034513		L	313028090		L	6034313		L	6034053		L	6034153		XL	ULNCS-4		L	3773	
	XL	6034304	10 boxes per carton/case	XL	6034514		XL	313028100	10 boxes per carton/case	XL	6034314		XL	6034054		XL	6034154		XXL	ULNCS-5		XL	3774	
	XXL	6034305																						

*See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

LATEX & VINYL RANGE



	MICROTOUCH® Plus			MICROTOUCH® NextStep®			MICROTOUCH® Plus Sterile			MICROTOUCH® E.P.®			MICROTOUCH® Elite®		
MATERIAL	Natural Rubber Latex			Natural Rubber Latex			Natural Rubber Latex			Natural Rubber Latex			Vinyl (Polyvinyl Chloride)		
KEY FEATURES & BENEFITS	<ul style="list-style-type: none">• Full textured surface for better grip• Inner polymer coating for easy donning• Excellent elasticity for a comfortable, natural fit• Superior tactile sensitivity and comfort			<ul style="list-style-type: none">• Full textured surface• Inner coating that hydrates and moisturizes skin• Excellent elasticity for a comfortable, natural fit			<ul style="list-style-type: none">• Superior strength and barrier protection• Inner polymer coating for easy donning• Fully textured surface for better grip• Packaged as a single sterile glove per pouch			<ul style="list-style-type: none">• Hand-specific for a comfortable, natural fit• Thicker cuff increases protection• Textured surface for sure handling of instruments• Extended 12" reinforced cuff			<ul style="list-style-type: none">• More comfortable and secure fit than traditional vinyl gloves• Smooth finish for enhanced tactile sensitivity• Protects from latex Type I allergy		
STERILE/NON-STERILE	Non-Sterile			Non-Sterile			Sterile			Non-Sterile			Non-Sterile		
GLOVE COLOR	Cream			Green			Cream			Cream			Cream		
CUFF LENGTH	Standard			Standard			Standard			Extended			Standard		
POWDER CONTENT	Powder-free			Powder-free			Powder-free			Powder-free			Powder-free		
EXTERNAL GLOVE SURFACE	Fully textured			Fully textured			Fully textured			Fully textured			Smooth		
FREEDOM FROM HOLES (AQL)	1.5			1.5			1.5			2.5			1.5		
CUFF THICKNESS	0.09 mm			0.1 mm			0.11 mm			0.5 mm			0.06 mm		
FINGER THICKNESS	0.16 mm			0.15 mm			0.17 mm			0.45 mm			0.1 mm		
PALM THICKNESS	0.13 mm			0.13 mm			0.15 mm			0.44 mm			0.11 mm		
GLOVE LENGTH	220 mm			240 mm			240 mm			300 mm			230 mm		
ALLERGY PREVENTION	N/A			N/A			N/A			N/A			Latex (Type 1)		
TESTED FOR USE WITH CHEMOTHERAPY DRUGS	No			No			No			No			No		
ANSELL TECHNOLOGIES	N/A			HYDRACARE™			N/A			N/A			N/A		
ORDERING INFORMATION	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING	SIZE	CODE	PACKAGING
	XS	6015300	150 gloves per box	XS	3200	100 gloves per box (size XS-L)	XS	N/A	1 glove / pouch 100 pouches per box 4 boxes per carton/case	XS	9201	1 pair / pouch 25 pouches per box 4 boxes per carton/case	XS	N/A	100 gloves per box
	S	6015301		S	3201	80 gloves per box (size XL)	S	6016001		S	9203		S	3091	
	M	6015302		M	3202		M	6016002		M	9205		M	3092	
	L	6015303	10 boxes per case/carton	L	3203	10 boxes per carton/case	L	6016003		L	9207		L	3093	10 boxes per case/carton
	XL	6015304		XL	3204		XL	N/A		XL	N/A		XL	3094	

CAUTION: This product contains natural rubber latex which may cause allergic reactions

SANDEL[®]

SAFETY SOLUTIONS

SHARPS SAFETY

Sharps Safety Solutions that help reduce the risk of sharps injuries and enhance compliance with recommendations from OSHA, AORN and ACS.

FLUID & TRIP MANAGEMENT

Simple, cost effective solutions for reducing the risk of slips, trips and falls for healthcare workers.

STAFF & PATIENT ERGONOMICS

Ergonomic solutions that are designed to reduce the risk of back and musculoskeletal injuries.

NEVER EVENT PROTECTION

Never Event Protection solutions help prevent the risk of costly surgical errors and aid in compliance with The Joint Commission.

PATIENT POSITIONERS

Patient positioners offer enhanced pressure relief and help reduce the risk of pressure injuries.

O.R. TURNOVER SOLUTIONS

Innovative, custom solutions to reduce turnover time, including proprietary antimicrobial linens which help reduce the risk of healthcare associated infections and improve patient outcomes.


SAFETY ESSENTIALS

A wide range of healthcare solutions designed to meet the unique requirements of healthcare facilities.

Change-A-Blade™ Safety Scalpel Handle

Designed for surgeon acceptance, Change-A-Blade™ is a disposable safety scalpel handle that provides the same weight, feel, visibility and puncture-ability as a traditional scalpel handle.

- Compatible with any manufacturer’s stainless steel or carbon blades
- Weighted to the feel of a conventional handle
- Allows for multiple blade changes during the same procedure
- Available in multiple handle sizes and accomodates several standard blade sizes
- Includes removable TIME OUT® Sleeve



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2200	Change-A-Blade™ Safety Scalpel Handle (#3)	Disposable, Weighted Safety Handle (#3) with TIME OUT Sleeve - Fits Any #10, #11 and #15 Blades (Sterile)	12/Bx 4 Bx/Cs	48
2200-L	Change-A-Blade™ Safety Scalpel Handle (#3-L)	Disposable, Weighted Safety Handle (#3-L) with TIME OUT Sleeve - Fits Any #10, #11 and #15 Blades (Sterile)	12/Bx 2 Bx/Cs	24
2204	Change-A-Blade™ Safety Scalpel Handle (#4)	Disposable, Weighted Safety Scalpel Handle (#4) with TIME OUT Sleeve - Fits Any #20, #22, #23 and #24 Blades (Sterile)	12/Bx 4 Bx/Cs	48

Available sterile and non-sterile




Weighted Safety Scalpels

SANDEL® Weighted Safety Scalpels mimic the weight and feel of a traditional scalpel handle. The scalpels feature a safety shield that retracts completely, providing similar visibility and puncture-ability to a traditional scalpel handle.

- Available in a variety of surgical blade and handle sizes
- Similar weight to a standard scalpel
- Color-coded safety shield for proper blade identification
- Suitable for right-handed and left-handed users
- Includes removable TIME OUT® Sleeve

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2210-N	Weighted Safety Scalpel	#3 Handle with #10 Blade, Orange Shield and TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2211-N	Weighted Safety Scalpel	#3 Handle with #11 Blade, Blue Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2215-N	Weighted Safety Scalpel	#3 Handle with #15 Blade, Yellow Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2220-N	Weighted Safety Scalpel	#4 Handle with #20 Blade, Orange Shield and TIME OUT Sleeve (Sterile)	12/Bx 4 Bx/Cs	48

Available sterile and non-sterile



384,000
sharp injuries are reported
in U.S. hospitals each year ¹

30%
of nurses suffer from a sharps
injury during their career ²

Non-Weighted Safety Scalpels

SANDEL® Non-Weighted Safety Scalpels feature a safety shield that retracts completely, providing similar visibility and puncture-ability to a traditional scalpel handle.

- Non-weighted – Ideal for ER, Cath Lab, Etc.
- Available in a variety of surgical blade sizes
- Suitable for right-handed and left-handed users
- Includes removable TIME OUT® Sleeve

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2230-N	Non-Weighted Safety Scalpel	#3 Handle with #10 Blade, Orange Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2231-N	Non-Weighted Safety Scalpel	#3 Handle with #11 Blade, Blue Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2235-N	Non-Weighted Safety Scalpel	#3 Handle with #15 Blade, Yellow Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96

Available sterile and non-sterile



Z-Friction Drape Neutral Zone

A non-magnetic instrument drape designated for use as a neutral zone for the hands-free transfer of sharps on the sterile field.

- Bright orange and highly recognizable
- Non-magnetic, non-slip design will help hold instruments where you place them



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1151	Z-Friction Drape Neutral Zone	Instrument Drape, Neutral Zone (10"x16") (Sterile)	40/Bx 2 Bx/Cs	80

Available sterile and non-sterile

Z-Tray & Stretch-A-Tray Hands-Free Transfer Trays

Disposable instrument trays designed to minimize the risk of sharps injuries through the hands-free transfer of sharps.

- Provides a designated area to safely place and retrieve sharp instruments
- Bright orange and highly recognizable
- Adhesive on bottom for the optional use as a non-movable neutral zone
- Stretch-A-Tray™ expands to accommodate a variety of instrument lengths



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1145	Z-Tray Hands-Free Transfer Tray	Transfer Tray with Adhesive Strip (Sterile)	12/Bx 2 Bx/Cs	24
1147	Stretch-A-Tray Hands-Free Transfer Tray	Extendable Transfer Tray with Adhesive Strip (Sterile)	10/Bx 2 Bx/Cs	20

Available sterile and non-sterile



1145




1147

SHARPS SAFETY

SharpsRink™ Neutral Zone

A non-magnetic, reusable neutral zone designed to support the hands-free transfer of sharps instruments.

- Bright orange and highly recognizable
- Non-magnetic, non-slip design will help hold instruments where you place them
- Ideal for facilities looking to reduce their environmental impact through the implementation of reusable products




ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1148-NS	SharpsRink™ Neutral Zone	Reusable Neutral Zone, 6" x 9" (Non-Sterile)	6/Bx 1 Bx/Cs	6

Disarm-It-All™


Counting, disarming and disposal boxes designed to safely remove scalpel blades, beaver blades and syringe needles.

- Allows users to safely remove scalpel blades, Beaver blades and syringe needles
- Two separate, removable halves with adhesive tape to provide stability on any surface
- Available with magnet or foam surface for counting needles and sutures
- Bright orange and highly recognizable




ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2243	Disarm-It-All™ (Magnet & Foam)	Foam Needle Counter (40 Count), Scalpel & Specialty Blade Remover, Syringe Needle Remover, Disposal Box	12/Bx 4 Bx/Cs	48
2244	Disarm-It-All™ (Magnet & Magnet)	Magnetic Needle Counter (30 Count), Scalpel & Specialty Blade Remover, Syringe Needle Remover, Disposal Box	12/Bx 4 Bx/Cs	48

Available sterile and non-sterile



2243



2244

Guard-It™ Syringe Capping Card


Provides a safe resting place for syringe needles and allows for a single handed technique for capping needles.

- Provides an alternative to the “scoop method”
- Compatible with any standard size syringe needle
- Bright orange and highly recognizable
- Disposable



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2280	Guard-It™	Syringe Capping Card (Sterile)	100/Bx 1 Bx/Cs	100

Available sterile and non-sterile



2280



FLUID & TRIP MANAGEMENT

Slips, trips & falls account for
25% of lost workday
injuries in
hospitals³

Trip-No-More™ Cord Cover

Bright orange cord cover designed to reduce the risk of trips and falls caused by cords and tubes on the floor.

- Serves as a “beacon” so healthcare personnel know where cords and tubes are located
- Adhesive strips keep it firmly in place
- Easy to apply, easy to remove
- Disposable, leaves no residue
- Bright orange and highly recognizable



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2302	Trip-No-More™ Cord Cover	Disposable Cord Cover (8" x 24") (Non-Sterile)	25/Bag 3 Bags/Cs	75
2303	Trip-No-More™ Cord Cover	Disposable Cord Cover (8" x 36") (Non-Sterile)	25/Bag 3 Bags/Cs	75
2309	Trip-No-More™ Cord Cover	Disposable Cord Cover (8" x 125' Roll) (Non-Sterile)	1 Roll/Cs	1
2312	Trip-No-More™ Cord Cover	Disposable Glow in the Dark Cord Cover (8" x 24") (Non-Sterile)	75/Bag 1 Bag/Cs	75
2313	Trip-No-More™ Cord Cover	Disposable Glow in the Dark Cord Cover (8" x 36") (Non-Sterile)	75/Bag 1 Bag/Cs	75
2350	Trip-No-More™ Cord Cover	Disposable Cord Cover (10" x 24") (Non-Sterile)	75/Bag 1 Bag/Cs	75
2351	Trip-No-More™ Cord Cover	Disposable Cord Cover (10" x 36") (Non-Sterile)	75/Bag 1 Bag/Cs	75
2360	Trip-No-More™ Cord Cover	Disposable Cord Cover (12" x 24") (Non-Sterile)	75/Bag 1 Bag/Cs	75
2361	Trip-No-More™ Cord Cover	Disposable Cord Cover (12" x 36") (Non-Sterile)	75/Bag 1 Bag/Cs	75
2369	Trip-No-More™ Cord Cover	Disposable Cord Cover (12" x 100' Roll) (Non-Sterile)	1 Roll/Cs	1



2302



2303



2309


DriFloor™ Absorbent Floor Pad

Absorbent floor pads designed to reduce the risk of slips and falls in the operating room or wherever fluids collect on the floor.


- Impervious, low-slip backing
- Super absorbent
- Eliminates the use of blankets and towels on the floor
- Helps reduce the risk of cross-contamination
- Helps improve turnover time
- Available in a variety of sizes




ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2419	DriFloor™ Absorbent Pad	Absorbent Floor Pad with Backing (16.5"x20") (Non-Sterile)	25/Bag 3 Bags/Cs	75
2440	DriFloor™ Absorbent Pad	Absorbent Floor Pad with Backing (33"x40") (Non-Sterile)	20/Bag 1 Bag/Cs	20
2446	DriFloor™ Absorbent Pad	Absorbent Floor Pad with Backing (35"x46") (Non-Sterile)	20/Bag 1 Bag/Cs	20
2425	DriFloor™ Absorbent Pad	Absorbent Floor Pad with Backing (28"x125") (Non-Sterile)	10/Bag 1 Bag/Cs	10
2449	DriFloor™ Absorbent Pad	Absorbent Floor Pad with Backing (33"x100') (Non-Sterile) - Roll	1 Roll/Cs	1




2446



2440



2419



2425

STAFF & PATIENTS ERGONOMICS

No-Lift Laws are enforced in
11 states⁶

80%
of nurses say they frequently
work with musculoskeletal pain⁴

The average cost of a patient handling-
related workers compensation claim is
\$15,600⁵

Z-Slider™ Patient Transfer Sheet

Disposable patient transfer sheet designed to prevent disabling back injuries to healthcare workers.

- Reduces the risk of physical strain to back, shoulders, neck and arms caused by moving patients
- No weight limit
- Single-patient use to reduce the risk of cross contamination
- Available in a variety of point-of-care dispensers for improved staff compliance
- Recyclable



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2101	Z-Slider™ - Single Pack	Patient Transfer Sheet, Individually Wrapped, 33" x 39" (Non-Sterile)	24/Bx 2 Bx/Cs	48
2151	Z-Slider™ XL	Patient Transfer Sheet, Individually Wrapped, 45" x 53" (Non-Sterile)	12/Bx 2 Bx/Cs	24
2102	Z-Slider™ - Roll	33" x 39" Patient Transfer Sheet in Roll Dispenser (Non-Sterile)	30/Bx 2 Bx/Cs	60
2103	Z-Slider™ Wall Dispenser (Large)	Patient Transfer Sheet Dispenser Refill for Large Point-of-Care Wall Holder (33" x 39" Patient Transfer Sheet) (Non-Sterile)	15/Bx 4 Bx/Cs	60
2104	Z-Slider™ Wall Dispenser (Small)	Patient Transfer Sheet Dispenser Refill for Small Point-of-Care Wall Holder (33" x 39" Patient Transfer Sheet) (Non-Sterile)	10/Bx 6 Bx/Cs	60
2101-H	Point-of-Care Clear Wall Holder	Point-of-Care Clear Wall Dispenser for use with Z-Slider™ Single Packs and Z-Slider™ XL Single Packs	4/Cs	4
2103-H	Point-of-Care Clear Wall Holder	Point-of-Care Clear Wall Mounting Holder for Use With 2103 Dispenser Box (12" x 20.5")	5/Cs	5
2104-H	Point-of-Care Clear Wall Holder	Point-of-Care Clear Wall Mounting Holder for the 2104 Dispenser Box (11.5" x 11")	5/Cs	5



2102



2103 & 2103-H



2101



2101-H



2104 & 2104-H

Z-TAP™ Patient Repositioning System

A system to support patient turning protocols and help prevent patient handling-related staff back injuries.


- Facilitates the Q2 turning protocol
- Highly breathable repositioning sheet and absorbent disposable underpad support skin integrity
- Low-friction repositioning sheet promotes safe patient handling



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
5101	Z-TAP™ Patient Repositioning Sheet	Patient Repositioning Sheet, 36" x 55"	20/case	20
5201	Z-SORB™ Disposable Underpad	Disposable Underpad, 33" x 51.2"	5/Bag, 10 bags/Cs	50
5301	Z-WEDGE™ Patient Positioning Wedges	30° Foam Wedges, 16" x 11" x 4.75", pair	1 Pr/Bag, 8 Bags/Cs	8 Pr/Cs



5101



5201



5301

ErgoPlus Anti-Fatigue Mats

Designed to reduce the stress and strain on muscles and joints caused from standing for extended procedures.

- Bounce-back memory foam restores energy and reduces strain on feet, legs and back
- Proven to reduce spinal compression and increase flexibility⁶
- Premium padding for additional comfort
- Available in various sizes
- Beveled edges help prevent tripping



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2317-P	ErgoPlus Anti-Fatigue Mat	Reusable Anti-Fatigue Mat (13" x 17")	5/Cs	5
2332-P	ErgoPlus Anti-Fatigue Mat	Reusable Anti-Fatigue Mat (20" x 32")	4/Cs	4
2348-P	ErgoPlus Anti-Fatigue Mat	Reusable Anti-Fatigue Mat (20" x 48")	4/Cs	4



2317-P



2332-P



2348-P

Ergo-Step Stools

Light-weight, stackable step stools which provide a safe working platform.

- Rubber top for added comfort
- Anti-skid feet keep stools firmly in place
- Can be stacked
- Interlocking feature creates a safe working platform
- Bright orange and highly recognizable
- NOTE: Top step stools (1170-T) should only be stacked on bottom step stools (1170-B). 1170-B can be stacked interchangeably as either a top or bottom stool



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1170	Ergo-Step Stool	2 Top Stools, 2 Bottom Stools (each 5.25"x18"x14")	4/Cs	4
1170-T	Ergo-Step Stool	4 Top Stools (each 5.25"x18"x14")	4/Cs	4
1170-B	Ergo-Step Stool	4 Bottom Stools (each 5.25"x18"x14")	4/Cs	4
1170-C	Ergo-Step Stool Cart	Step Stool Cart - Holds a Maximum of 6 Stools	1/Cs	1



1170-T



1170



1170-B



1170-C

Sit-Stand Stool

A stool designed to help relieve pressure and fatigue during long procedures.

- Ideal for the O.R. and nurses stations
- Height 20" (min.) to 30" (max.)
- Weight capacity – 330 lbs



2300

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2300	Sit-Stand Stool	Sit-Stand Stool (Non-Sterile)	1/Cs	1

The annual cost of measurable medical errors in the U.S. has been estimated at **\$17.1 billion**⁷

The Joint Commission requires compliance with the National Patient Safety Goals for medication safety and Universal Protocol for preventing wrong-site surgery⁸

Correct Medication Labeling System

The SANDEL® Correct Medication Labeling System provides an easy way to comply with AORN guidelines and Joint Commission standards for medication safety. Choose from a wide variety of components including preprinted labels, color-coding flags, TIME OUT® markers and the SANDEL® Specimen Zone.



MEDICATION LABELS

- Preprinted labels
- Moisture and smear resistant
- Will not tear or fall off when submerged in fluid

COLORED ID FLAGS WITH SYRINGE ID STRIP

- Proprietary flag system increases visibility of labels
- Flags adhere to any surface yet leave no residue
- Includes matching peel-off syringe ID strip for easy identification

MARKERS

- A variety of skin & permanent markers to select from
- “No-roll” triangular barrel

SPECIMEN ZONE

- Assists in organizing and labeling specimens on the sterile field
- Non-woven, perforated sheet
- Pre-numbered
- Fluid-proof backing helps keep specimen moist while on the sterile field


OPERATING ROOM

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
3302	Correct Medication Labeling System - O.R.	(4) Flags, (1) 4-in-1 Marker, (2) Sheets of 24 Preprinted Labels Designed for the O.R., 2 Med Cups, (1) Specimen Zone (4 count) (Sterile)	25/Bx 2 Bx/Cs	50
3303-P	Correct Medication Labeling System - O.R.	(4) Flags, (1) 4-in-1 Marker, (2) Sheets of 24 Preprinted Labels Designed for the O.R., (1) Specimen Zone (4 count) (Sterile)	50/Bx 2 Bx/Cs	100
3304	Correct Medication Labeling System - O.R.	(4) Flags, (2) Sheets of 24 Preprinted Labels Designed for the O.R., (1) Specimen Zone (4 count) (Sterile)	50/Bx 1 Bx/Cs	50
3310-P	Correct Medication Labeling System - O.R.	(2) Flags, 4-in-1 Marker, (2) Sheets of 24 Preprinted Labels Designed for the O.R., (1) Specimen Zone (4 count) (Sterile)	50/Bx 2 Bx/Cs	100
3311	Correct Medication Labeling System - O.R.	(2) Flags, (1) Petite Permanent Marker, (2) Sheets of 16 Preprinted Labels Designed for the O.R. (Sterile)	100/Bx 1 Bx/Cs	100
3312	Correct Medication Labeling System - O.R.	(2) Flags, (2) Sheets of 24 Preprinted Labels Designed for the O.R., (1) Specimen Zone (4 count) (Sterile)	50/Bx 1 Bx/Cs	50
3330-P	Correct Medication Labeling System - O.R.	(2) Sheets of 24 Preprinted Labels Designed for the O.R. (Sterile)	100/Bx 2 Bx/Cs	200
3331-P	Correct Medication Labeling System - O.R.	(2) Sheets of 16 Preprinted Labels Designed for the O.R. (Sterile)	100/Bx 2 Bx/Cs	200
3333	Correct Medication Labeling System - O.R.	(2) Sheets of 16 Preprinted Labels Designed for the O.R., (1) Petite Permanent Marker (Sterile)	100/Bx 1 Bx/Cs	100
3334-P	Correct Medication Labeling System - O.R.	(2) Sheets of 24 Preprinted Labels Designed for the O.R., (1) 4-in-1 Marker, (1) Specimen Zone (Sterile)	50/Bx 2 Bx/Cs	100

Available sterile and non-sterile



Skin & Permanent Markers



Specimen Zone



Medication Cup

Antibiotic Strength: _____ Expiration: _____ Date / Time: _____	Saline - Injectable Strength: 0.9% Expiration: _____ Date / Time: _____
Bacitracin Strength: _____ Expiration: _____ Date / Time: _____	Saline - Irrigation Strength: 0.9% Expiration: _____ Date / Time: _____
Epinephrine Strength: _____ Expiration: _____ Date / Time: _____	Water Strength: _____ Expiration: _____ Date / Time: _____
Heparin Strength: _____ Units Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Lidocaine Strength: _____ % Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Lidocaine & Epinephrine Strength: _____ % Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Marcaine Strength: _____ % Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Marcaine & Epinephrine Strength: _____ % Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____

OPERATING ROOM - 16 Sandel Medical Industries, LLC PL No 120713 Rev. B
Use Permanent Ink Only - Do Not Use Gelatin Viallet Ink

Operating Room (16)

Afrin Strength: _____ Expiration: _____ Date / Time: _____	Marcaine Strength: _____ % Expiration: _____ Date / Time: _____
Antibiotic - 0.9% Saline Irrigation Strength: _____ Expiration: _____ Date / Time: _____	Marcaine & Epinephrine Strength: _____ % Expiration: _____ Date / Time: _____
Bacitracin Strength: _____ Units Expiration: _____ Date / Time: _____	Papaverine Strength: _____ Expiration: _____ Date / Time: _____
Cocaine Strength: _____ % Expiration: _____ Date / Time: _____	Saline - Injectable Strength: 0.9% Expiration: _____ Date / Time: _____
Contrast Strength: _____ Expiration: _____ Date / Time: _____	Saline - Irrigation Strength: _____ Expiration: _____ Date / Time: _____
Epinephrine Strength: _____ Expiration: _____ Date / Time: _____	Thrombin Strength: _____ Units Expiration: _____ Date / Time: _____
Heparin Strength: _____ Units Expiration: _____ Date / Time: _____	Water Strength: _____ Expiration: _____ Date / Time: _____
Heparin Strength: _____ Units Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Kefzol Strength: _____ Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Lidocaine Strength: _____ % Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Lidocaine & Epinephrine Strength: _____ % Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____
Med. Strength: _____ Expiration: _____ Date / Time: _____	Med. Strength: _____ Expiration: _____ Date / Time: _____

OPERATING ROOM Sandel Medical Industries, LLC PL No 120714 Rev. A
Use Permanent Ink Only - Do Not Use Gelatin Viallet Ink

Operating Room (24)

Correct Medication Labeling System

SURGERY CENTER

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
3300	Correct Medication Labeling System - Surgery Centers	(2) Flags, (1) 4-in-1 Marker, (1) Sheet of 16 Preprinted Labels Designed for Surgery Centers, (1) Specimen Zone (4 count) (Sterile)	50/Bx 1 Bx/Cs	50

Available sterile and non-sterile

GENERAL USE

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
3300-1	Correct Medication Labeling System - General Use	(2) Flags, (1) 4-in-1 Marker, (1) Sheet of 16 Preprinted Labels with Prompts for Medication Name and Strength. Designed for General Use (Sterile)	50/Bx 1 Bx/Cs	50

Available sterile and non-sterile

OTHER

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
3301	Correct Medication Labeling System - Cath Labs	(4) Flags, (1) 4-in-1 Marker, (2) Sheets of 16 Preprinted Labels Designed for Cath Labs (Sterile)	50/Bx 1 Bx/Cs	50
3305	Correct Medication Labeling System – Flags Only	(4) Flags (1 ea. Orange, Red, Yellow, Blue) (Sterile)	50/Bx 1 Bx/Cs	50
3350	Correct Medication Labeling System - Interventional Radiology	(2) Flags, (1) 4-in-1 Marker, (2) Sheets of 16 Preprinted Labels Designed for Interventional Radiology (Sterile)	50/Bx 1 Bx/Cs	50
3360	Correct Medication Labeling System - Ophthalmic	(1) Flag, (1) Permanent Marker, (1) Sheet of 16 Preprinted Labels Designed for Ophthalmic (Sterile)	50/Bx 1 Bx/Cs	50
3370	Correct Medication Labeling System - Pain Management	(2) Flags, (1) 4-in-1 Marker, (2) Sheets of 16 Preprinted Labels Designed for Pain Management (Sterile)	50/Bx 1 Bx/Cs	50
3390	Correct Medication Labeling System - Labor & Delivery	(2) Flags, (1) 4-in-1 Marker, (1) Sheet of 16 Preprinted Labels Designed for Labor & Delivery (Sterile)	50/Bx 1 Bx/Cs	50
3334-P	Correct Medication Labeling System - O.R.	(2) Sheets of 24 Preprinted Labels Designed for the O.R., (1) 4-in-1 Marker, (1) Specimen Zone (Sterile)	50/Bx 2 Bx/Cs	100

Available sterile and non-sterile

CUSTOM LABELS

To receive a quote for a custom medication labeling system, or custom labels contact SANDEL Customer Service at 866-764-3327 or sandel-cs@ansell.com.

Correct Medication Labeling System

Antibiotic 0.9% Saline Irrigation Name and Strength Expiration Date / Time	Antibiotic 0.9% Saline Irrigation Name and Strength Expiration Date / Time
Lidocaine Strength: _____ % Expiration Date / Time	Lidocaine Strength: _____ % Expiration Date / Time
Lidocaine E Epinephrine Strength: _____ % Expiration Date / Time	Lidocaine E Epinephrine Strength: _____ % Expiration Date / Time
Marcaine Strength: _____ % Expiration Date / Time	Marcaine Strength: _____ % Expiration Date / Time
Marcaine E Epinephrine Strength: _____ % Expiration Date / Time	Marcaine E Epinephrine Strength: _____ % Expiration Date / Time
Saline - Irrigation Strength: 0.9% Expiration Date / Time	Saline - Irrigation Strength: 0.9% Expiration Date / Time
Med. Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Med. Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time

SLURBURY CENTER Sandel Medical Industries, LLC PI No. 120030 Rev. A
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

Surgery Centers

Antibiotic 0.9% Saline Irrigation Name and Strength Expiration Date / Time	Antibiotic 0.9% Saline Irrigation Name and Strength Expiration Date / Time
Lidocaine Strength: _____ % Expiration Date / Time	Lidocaine Strength: _____ % Expiration Date / Time
Lidocaine E Epinephrine Strength: _____ % Expiration Date / Time	Lidocaine E Epinephrine Strength: _____ % Expiration Date / Time
Marcaine Strength: _____ % Expiration Date / Time	Marcaine Strength: _____ % Expiration Date / Time
Marcaine E Epinephrine Strength: _____ % Expiration Date / Time	Marcaine E Epinephrine Strength: _____ % Expiration Date / Time
Saline - Irrigation Strength: 0.9% Expiration Date / Time	Saline - Irrigation Strength: 0.9% Expiration Date / Time
Med. Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Med. Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time

GENERIC Sandel Medical Industries, LLC Patent Pending 120044 Rev. A
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

General Use

Adenosine Strength: _____ Expiration Date / Time	Lidocaine Strength: _____ % Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Nitroglycerin Strength: _____ Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Papaverine Strength: _____ Expiration Date / Time
Heparin Strength: _____ Units Expiration Date / Time	Verapamil Strength: _____ Expiration Date / Time
Heparin 0.9% Saline Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Heparin 0.9% Saline Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Heparin 0.9% Saline Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time

CATH. LAB. Sandel Medical Industries, LLC PI No. 120043 Rev. A
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

Cath Labs

Contrast Strength: _____ Expiration Date / Time	Heparin Strength: _____ Units Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Heparin 0.9% Saline Strength: _____ Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Lidocaine Strength: _____ % Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Lidocaine Strength: _____ % Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Nitroglycerin Strength: _____ Expiration Date / Time
Heparin Strength: _____ Units Expiration Date / Time	Sodium Bicarbonate Strength: _____ Expiration Date / Time
Heparin 0.9% Saline Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Heparin 0.9% Saline Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time

INTERVENTIONAL RADIOLOGY Sandel Medical Industries, LLC PI No. 120042 Rev. A
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

Interventional Radiology

Air Strength: _____ Expiration Date / Time	Lidocaine / Xylocaine E Epinephrine Strength: _____ % Expiration Date / Time
Antibiotic Strength: _____ Expiration Date / Time	Marcaine Strength: _____ % Expiration Date / Time
B S S Strength: _____ Expiration Date / Time	Miochol Strength: _____ Expiration Date / Time
Cornea Coat Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Epinephrine Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Gentamycin Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Indocyanine Green Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Lidocaine / Xylocaine Strength: _____ % Expiration Date / Time	Med. Strength: _____ Expiration Date / Time

OPHTHALMIC Sandel Medical Industries, LLC PI No. 120481 Rev. A
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

Ophthalmic Surgery

Ancef Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Contrast Strength: _____ Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
DepoMedrol Strength: _____ mg Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Kenalog Strength: _____ mg Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Lidocaine Strength: _____ % Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
% Lidocaine E % Sodium Bicarbonate Expiration Date / Time	
Marcaine Strength: _____ % Expiration Date / Time	
Marcaine E Epinephrine Strength: _____ % Expiration Date / Time	

PAIN MANAGEMENT Sandel Medical Industries, LLC PI No. 120478 Rev. B
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

Pain Management

Hemabate Strength: _____ Expiration Date / Time	Water Irrigation Strength: _____ Expiration Date / Time
Lidocaine Strength: _____ % Expiration Date / Time	Water Irrigation Strength: _____ Expiration Date / Time
Lidocaine Strength: _____ % Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Nesacaine Strength: _____ % Expiration Date / Time	Med. Strength: _____ Expiration Date / Time
Nesacaine Strength: _____ % Expiration Date / Time	
Normal Saline Strength: 0.9% Irrigation Expiration Date / Time	Betadine
Normal Saline Strength: 0.9% Irrigation Expiration Date / Time	Betadine
Pitocin Strength: _____ Units Expiration Date / Time	Hibiclens

LABOR & DELIVERY Sandel Medical Industries, LLC PI No. 120592 Rev. A
Use Permanent Ink Only - Do Not Use Geriatric Violet Ink

Labor & Delivery

PREP-RESIST™ Skin Markers

SANDEL PREP-RESIST™ skin markers are resistant to chlorhexidine gluconate and isopropyl alcohol skin prep solutions. Specially formulated to show up on varying colors of skin, PREP-RESIST markers offer good visibility after prep but largely fade after a day so patients aren't burdened with lingering surgical site markings.

- Resistant to chlorhexidine gluconate and isopropyl alcohol skin prep solutions
- Gentian violet-free
- Available in standard length, petite, and 4-in-1 dual-tip/dual-ink marker with chiseled tip
- Triangular “no-roll” barrel
- Bright orange and highly recognizable
- 4-in-1 marker comes with TIME OUT® sleeve



Post Prep: PREP-RESIST Marker vs. Traditional Skin Marker

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1015-NNS	PREP-RESIST™ Petite Skin Marker	Petite Prep-Resistant “No Roll” Skin Marker (Non-Sterile)	50/Bag 10 Bags/Cs	500
1029	PREP-RESIST™ Skin Marker	Prep-Resistant “No Roll” Skin Marker (Sterile)	25/Bag 4 Bx/Cs	100
1029-NNS	PREP-RESIST™ Skin Marker	Prep-Resistant, “No Roll” Skin Marker (Non-Sterile)	200/Bag 2 Bags/Cs	400
1061	PREP-RESIST™ 4-in-1 Marker	Dual Tip/Dual Ink Marker (Permanent and Prep-Resistant Skin Ink) with Chiseled Tip, TIME OUT Sleeve (Sterile)	25/Bx 4 Bx/Cs	100
1061-NNS	PREP-RESIST™ 4-in-1 Marker	Dual Tip/Dual Ink Marker (Permanent and Prep-Resistant Skin Ink) with Chiseled Tip, TIME OUT Sleeve (Non-Sterile)	100/Bag 4 Bags/Cs	400



TIME OUT® Skin Markers & Permanent Markers

TIME OUT® Skin Markers are designed to mark the correct surgical site on a patient’s skin prior to a procedure and in compliance with The Joint Commission’s Universal Protocol.

- Gentian violet-free
- Triangular barrel prevents rolling
- Available in single and dual-tip
- SANDEL 4-in-1 Marker features a chiseled tip to allow the user to make lines with multiple thickness options
- Permanent marker for medication labeling and other purposes
- Bright orange and highly recognizable
- Disposable




ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1000-NNS	Permanent Marker	“No-Roll” Permanent Marker, Black Ink (Non-Sterile)	100/Bag 4 Bags/Cs	400
1002-NNS	Petite Permanent Marker	Petite “No-Roll” Permanent Marker, Black Ink (4” in Length) (Non-Sterile)	50/Bag 10 Bags/Cs	500
1011-NNS	Petite Skin Marker	Petite “No Roll” Skin Marker (4” in Length) (Non-Sterile)	50/Bag 10 Bags/Cs	500
1011-NS1	Petite Skin Marker	Petite “No Roll” Skin Marker (4” in Length) (Non-Sterile)	500/Bx 2 Bx/Cs	1000
1019	Skin Marker	“No Roll” Skin Marker (Sterile)	25/Bx 4 Bx/Cs	100
1019-R	Skin Marker w/Flexible Ruler	“No Roll” Skin Marker with Flexible Ruler (Sterile)	25/Bx 4 Bx/Cs	100
1019-NNS	Skin Marker	“No Roll” Skin Marker (Non-Sterile)	200/Bag 2 Bags/Cs	400
1031	2-in-1 Marker	“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), TIME OUT Sleeve (Sterile)	25/Bx 4 Bx/Cs	100
1031-LR	2-in-1 Marker w/Labels & Flexible Ruler	“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT Sleeve & Flexible Ruler (Sterile)	25/Bx 4 Bx/Cs	100
1031-LRNS	2-in-1 Marker w/Labels & Flexible Ruler	“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT Sleeve & Flexible Ruler (Non-Sterile)	100/Bag 2 Bags/Cs	200
1031-NNS	2-in-1 Marker	“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), TIME OUT Sleeve (Non-Sterile)	200/Bag, 2 Bags/Cs	200
1041	4-in-1 Marker	“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), Skin Marker with Chiseled Tip, TIME OUT Sleeve (Sterile)	25/Bx 4 Bx/Cs	100
1041-R	4-in-1 Marker with Ruler	“No Roll” Dual Tip/Dual Ink (Permanent and Skin Ink) Skin Marker with Chiseled Tip, TIME OUT Sleeve and Flexible Ruler (Sterile)	25/Bx, 4 Bx/Cs	100
1041-L	4-in-1 Marker w/Labels	“No Roll” Dual Tip/Dual Tip (Permanent & Skin Ink) Marker with 8 labels & TIME OUT Sleeve (Sterile)	25/Bx, 4 Bx/Cs	100
1041-LR	4-in-1 Marker w/Labels & Flexible Ruler	“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT Sleeve & Flexible Ruler (Sterile)	25/Bx 4 Bx/Cs	100
1041-LRNS	4-in-1 Marker	“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT Sleeve & Flexible Ruler (Non-Sterile)	100/Bag 2 Bags/Cs	200
1041-NNS	4-in-1 Marker	“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), Skin Marker with Chiseled Tip, TIME OUT Sleeve (Non-Sterile)	100/Bag 4 Bags/Cs	400
1041-NS1	4-in-1 Marker	“No-Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink) Skin Marker with Chiseled Tip (Sterile)	100/bag, 4 bags/Cs	400

Correct Site Stickers and Tattoos

Correct Site stickers and tattoos offer an easy way to mark the surgical site and comply with Joint Commission and AORN recommendations for correct site marking.

- Ideal for use on dark-skinned patients
- Convenient size
- Tattoos are legible for up to 7 days
- Latex-free



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
3176	Correct Site Sticker	Small (1" diameter) Green and White Circle Sticker (Non-Sterile)	100/Bag, 1 Bag/Cs	100
3186	Correct Site Tattoo	Small (1" diameter) Green & White Circle Tattoo (Non-Sterile)	50/Bag, 5 Bags/Cs	250

TIME OUT® Sleeve

The TIME OUT® Sleeve is designed to be placed over a scalpel as a reminder to comply with The Joint Commission’s Universal Protocol.

- Bright orange and highly recognizable
- Ruler on side panel
- Included with the SANDEL® Change-A-Blade™ Safety Scalpel Handle, Weighted and Non-Weighted Safety Scalpels, 4-in-1 and 2-in-1 Markers



911

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
911	TIME OUT® Sleeve	Disposable TIME OUT Sleeve (Sterile)	300/Bx 1 Bx/Cs	300

Available sterile and non-sterile

TIME OUT® Beacon

TIME OUT® Beacons serve as a prominent reminder for the surgical team to comply with The Joint Commission’s Universal Protocol.

- Can be placed over the surgical site, mayo stand or back table
- Bright orange and highly recognizable
- Available in cotton or a non-woven biodegradable material



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
811	TIME OUT® Beacon - Cotton	Disposable TIME OUT Reminder (Cotton) (Sterile)	24/Bx 4 Bx/Cs	96
811-E	TIME OUT® Beacon - Non-Woven	Disposable TIME OUT Reminder (Non-Woven) (Sterile)	100/Bx 2 Bx/Cs	200

Available sterile and non-sterile



811



811-E

COUNT IN PROGRESS Beacon

The Count in Progress Beacon is designed to help reduce the risk of retained surgical items in the O.R. by eliminating distractions during the counting process.

- Notifies the surgical team that a count is about to begin
- Helps reduce distractions and unnecessary activity while the count is in progress
- Bright orange and highly recognizable



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
812-C	Count In Progress Beacon (Sterile)	Count In Progress Beacon - Cotton (Sterile)	24/Bx 4Bx/Cs	96

Available sterile and non-sterile


Specimen Handling System

The Specimen Handling System aids compliance with AORN recommendations for specimen management, including labeling, identification, and protection during transport. With Specimen Handling System products, specimens are easily recognized to help prevent accidental disposal or other mismanagement.

- Leak-proof, puncture resistant specimen cups
- Moisture-proof orange labels
- Bright orange, biohazard labeled bags to protect workers from infectious specimens
- Easy-to-handle heart shaped specimen bowl
- Specimen Zone aids identification of tissue samples
- All items sold separately



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
3450	Correct Specimen Zone	Orange Specimen Zone (4 Count) (Sterile)	50/Bx 2 Bx/Cs	100
3455	Orange Specimen Labels	Orange Specimen Labels (2"x4") on roll (Non-Sterile)	500/Roll 1 Roll/Cs	500
3460	Orange Specimen Wrap	Orange Specimen Wrap (23 1/2" x 16 1/2") (Sterile)	24/Bx 4 Bx/Cs	96
3460-NS	Orange Specimen Wrap	Orange Specimen Wrap (23 1/2" x 16 1/2") (Non-Sterile)	100/Bag 1 Bag/Cs	100
3470	Orange 120cc Cup with Lid	120cc Cup with Orange Lid (Sterile)	25/Bx 2 Bx/Cs	50
3475	Orange Specimen Bowl	Shaped Orange Specimen Bowl (Sterile)	20/Bag 1 Bags/Cs	20
3475-NS	Heart Shaped Specimen Bowl	Heart Shaped Specimen Bowl (Non-Sterile)	100/Bag, 1 Bag/Cs	100
3480	Orange Specimen Bag	Orange Specimen Bag Large (33 3/8" x 39 1/4") (Non-Sterile)	25/Bag 2 Bags/Cs	50



3450



3455



3460



3470



3475



3480

SANDEL® Safety Kits

SANDEL® Safety Kits offer a wide variety of solutions designed to fit the unique needs of your facility. Choose from an existing safety kit or create your own.



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
K-911	TIME OUT® Safety Kit	(1) Friction Drape (10 x 16) (1) Stretch-A-Tray – Hands-Free Transfer (1) Specimen Zone (4 count) (1) Correct Medication Labeling System: 4 color flags, 2 sheets of 24 preprinted labels, 1 ea. 4-in-1™ Marker (Sterile)	12/Bx 2 Bx/Cs	24
K-911-10	TIME OUT® Safety Kit	(1) #10 Weighted Safety Scalpel with TIME OUT Sleeve (1) Z-Friction Drape (10" x 16") (1) Z-Tray – Hands-Free Transfer (1) Specimen Zone (4 Count) (1) Correct Medication Labeling System: 4 color flags, 2 sheets of 24 preprinted labels, 1 ea. 4-in-1™ Marker (Sterile)	12/Bx 2 Bx/Cs	24
K-911-15	TIME OUT® Safety Kit	(1) #15 Weighted Safety Scalpel with TIME OUT Sleeve (1) Z-Friction Drape (10" x 16") (1) Z-Tray – Hands-Free Transfer (1) Specimen Zone (4 Count) (1) Correct Medication Labeling System: 4 color flags, 2 sheets of 24 preprinted labels, (1) 4-in-1 Marker (Sterile)	12/Bx 2 Bx/Cs	24
2271-B	Change-A-Blade® Safety Kit	(1) Change-A-Blade, (1) Disarm-It-All with Magnet & Foam (2) Sheets of 16 O.R. labels, (1) 4-in-1 Marker (4) Color Flags, TIME OUT Beacon - Woven (Sterile)	10/Bx 4 Bx/Cs	40


The cost to treat a single pressure ulcer can reach up to **\$150,000.**⁹

ProForm™ Patient Positioners

Help reduce the risk of pressure injuries in your O.R. ProForm Patient Positioners offer enhanced pressure relief compared to traditional foam positioning products.¹⁰ Each ProForm patient positioner reduces pressure by 20% - 50% compared to a traditional foam patient positioner.¹¹

- 1. Anatomical shape
 - Improved weight distribution and pressure relief
- 2. Supine head positioners and ulnar nerve protectors are made of memory foam
- 3. Disposable
 - Helps reduce the risk of cross-contamination and improve O.R. turnover time
- 4. Bright orange
 - Highly recognizable in the operating room

ITEM #	PRODUCT NAME	DESCRIPTION	CASE QTY.
4100	ProForm Patient Positioner	Supine Head Positioner - Contoured - Memory Foam	24/Cs
4120	ProForm Patient Positioner	Supine Head Positioner - Memory Foam	36/Cs
4150	ProForm Patient Positioner	Prone Head Positioner - Resilient Foam	24/Cs
4200	ProForm Patient Positioner	Ulnar Nerve Protectors - Memory Foam	36 pr/Cs
4300	ProForm Patient Positioner	Arm Cradles - Laminectomy - Resilient Foam	12 pr/Cs
4400	ProForm Patient Positioner	Heel Positioners - Resilient Foam	36 pr/Cs
4500	ProForm Patient Positioner	45" O.R. Body Wedge , 18" x 8" x 8" - Resilient Foam	8/Cs



Healthcare associated infections (HAIs) occur nearly **2 million times** each year.¹²

HAIs cost **\$36 to \$45 billion** annually in the US.¹³

STAT-PAC™ Turnover Packs

Standard or custom packs are available with a variety of components to meet the needs of any department or specialty including the operating room, labor and delivery, cystoscopy, catheterization lab, radiology, endoscopy, and more.

- Improve consistency, efficiency and turnover time
- Enhance infection control
- Decrease lint to help reduce the risk of cross-contamination
- Help maintain patient skin integrity
- Eliminate costs related to reusable textiles
- Comply with guidelines from AORN
- Orange color is a visible indicator of the antimicrobial and provides separation from blue sterile drapes
- Extensive offering of turnover pack components to choose from – additional options available upon request



ITEM #	PRODUCT NAME	DESCRIPTION
SP-02	STAT-BLOC Linen Set (ORANGE)	Antimicrobial, Absorbent Table Sheet, Lift Sheet, Armboard Covers
SP-03	STAT-BLOC Linen Set (ORANGE)	Antimicrobial, Absorbent Table Sheet, Lift Sheet, Armboard Covers, Headrest Cover

Linen sets can be included in custom turnover packs with other components such as bags, patient positioning straps, mops and wipes.



SP-02



SP-03

STAT-BLOC™ Table / Gurney Sheets

Disposable, antimicrobial, low-linting table sheets provide three layers of protection. Suitable for operating room tables, cystoscopy tables, gurneys, stretchers and beds.

- 1. Quilted, breathable top layer allows for fluid absorption and minimizes wrinkling to help prevent skin breakdown
- 2. Strong concentration of Super Absorbent Polymers (SAP) for maximum fluid control to wick fluid away from the patient's skin (can absorb up to 5 liters of fluid*)
- 3. Impervious, antimicrobial backing prevents bacteria from colonizing on the sheet
 - Kills 99% of COVID-19 which comes into contact with the antimicrobial backing*
 - Antimicrobial is proven 99.9% effective against MRSA, CRE and E.Coli*



ITEM #	PRODUCT NAME	DESCRIPTION
AS-92-QZ	STAT-BLOC™ Table Sheet	General Surgery - Disposable, Antimicrobial, Absorbent, Quilted Table Sheet
AS-100-QZ	STAT-BLOC™ Table Sheet	Cath Lab / Interventional Radiology - Disposable, Antimicrobial, Absorbent, Quilted, XL Table Sheet
AS-92-QZF	STAT-BLOC™ Table Sheet	Transport / Stretcher Sheet - Disposable, Antimicrobial, Absorbent, Quilted, Fitted Sheet
AS-92-QXLF	STAT-BLOC™ Table Sheet	Endoscopy / Ophthalmic - Disposable, Antimicrobial, Absorbent, Quilted, Fitted Table Sheet with XL Absorbent Area
AS-6030-QZ	STAT-BLOC™ Table Sheet	Cystoscopy - Disposable, Antimicrobial, Absorbent, Quilted, Top and Bottom Table Sheets
CF-AS9-LS-QZP-C	STAT-BLOC™ Table Sheet	Pediatrics - Disposable, Antimicrobial, Absorbent, Quilted Table Sheet, Orange with Sunshine Pattern (Includes standard Draw Sheet)
WGP-01	STAT-BLOC™ Table Sheet / Cover	Pediatrics - Disposable, Antimicrobial Cover for Gurney or Wagon, Orange with Sunshine Pattern

* Data on file.

SANDEL® Draw Sheets

Disposable, non-woven, low linting draw sheets for transferring and repositioning patients.

- Draw sheets are available combo-folded with table sheets or separately
- Standard draw sheet is capable of supporting up to 1,000 pounds of weight in both wet and dry conditions*



ITEM #	PRODUCT NAME	DESCRIPTION
LFS-66	Draw Sheet	Disposable, Standard Draw Sheet, White
LFS-4080	Draw Sheet	Disposable, XL Draw Sheet, White - Bariatrics / Cath Lab / Interventional Radiology
LFS-66-BF	Draw Sheet	Disposable, Butterfly Draw Sheet, White

STAT-BLOC™ Head Rest, Pillow and Stirrup Covers

Disposable, Antimicrobial Stirrup Covers for use on pillows, headrests and stirrups to help reduce the risk of cross-contamination.

- Impervious, antimicrobial layer prevents bacteria from colonizing on the cover
- Antimicrobial is proven 99.9% effective against MRSA, CRE and E.Coli*




ITEM #	PRODUCT NAME	DESCRIPTION
DPC-10	STAT-BLOC™ Pillow Cover	Disposable, Antimicrobial Pillow Cover
HR-10	STAT-BLOC™ Headrest Cover	Disposable, Antimicrobial Headrest Cover
SC-01	STAT-BLOC™ Stirrup Cover	Disposable, Antimicrobial Stirrup Cover



DPC-10



HR-10



SC-01

* Data on file.

STAT-BLOC™ Arm Board Covers

Disposable, antimicrobial covers for use on arm boards to help reduce the risk of cross-contamination.

- Impervious, antimicrobial layer prevents bacteria from colonizing on the cover
- Antimicrobial is proven 99.9% effective against MRSA, CRE and E.Coli*
- Dual straps for added safety
- Straps secured to cover to prevent straps from falling to the floor



ITEM #	PRODUCT NAME	DESCRIPTION
ABC-10	STAT-BLOC™ Arm Board Cover	Disposable, Antimicrobial Arm Board Cover with 2 Straps
ABC-XL	STAT-BLOC™ Arm Board Cover	Disposable, Antimicrobial, XL Arm Board Cover with 2 Straps



ABC-10



ABC-XL

STAT-STRAP™ Patient Positioning Straps

Disposable straps used to safely and comfortably position patients on the operating room table, cystoscopy tables, gurneys, stretchers and beds.

- Disposable to minimize the risk of healthcare-associated infections
- Fit around all areas of the body, including the thighs and legs



ITEM #	PRODUCT NAME	DESCRIPTION
TS-SM	STAT-STRAP™	Disposable Patient Positioning Strap, Standard, Velcro Closure, 2 piece
TS-10	STAT-STRAP™	Disposable Patient Positioning Strap, Long, Velcro Closure, 2 piece
TS-SS	STAT-STRAP™	Disposable Patient Positioning Strap, Soft, Velcro Closure, 2 piece
TS-PD	STAT-STRAP™	Disposable Patient Positioning Strap, Padded, Velcro Closure, 2 piece



TS-SM / TS-10



TS-SS



TS-PD

STAT-MOP™

A specially designed mop for use in the operating room to provide higher absorbency and decrease lint.

- Disposable to help prevent cross-contamination
- More absorbent than traditional cotton, rayon and microfiber mops¹⁴
- Low-linting
- Made from 100% post-industrial material*
- Meets AORN guidelines for using a clean mop between each case¹⁵

ITEM #	PRODUCT NAME	DESCRIPTION
STM-MD	STAT-MOP™	Standard Mop Head
STM-MD-RF	STAT-MOP™	Standard Mop Head with Fan Tail



STM-MD



STM-MD-RF

Rayon Mops

Traditional disposable mops for cleaning the floors during room turnover.

- Disposable to help prevent cross-contamination
- Meets AORN guidelines for using a clean mop between each case¹⁵

ITEM #	PRODUCT NAME	DESCRIPTION
M16-R	Mop Head	Rayon Mop Head 12 oz. (#16)
M16-RF	Mop Head	Rayon Mop Head 12 oz. (#16) with Fan Tail
M14OZ-R	Mop Head	Rayon Mop Head 14 oz.
M14OZ-RF	Mop Head	Rayon Mop Head 14 oz. with Fan Tail
M16OZ-R	Mop Head	Rayon Mop Head 16 oz.
M16OZ-RF	Mop Head	Rayon Mop Head 16 oz. with Fan Tail



M16-R



M16-RF

* Data on file.

Microfiber Mops

Disposable mops designed to attract and trap dust, dirt and bacteria from the surface for superior cleaning while using less chemicals and water than traditional mops.

- Disposable to help prevent cross-contamination
- Highly absorbent
- Non-linting
- Meets AORN guidelines for using a clean mop between each case¹⁵

ITEM #	PRODUCT NAME	DESCRIPTION
MMB-T479RM-18	Mop Head	Microfiber Flat Mop, Blue
MMW-616-P	Mop Head	Microfiber & Cotton Flat Mop, White



MMB-T479RM-18



MMW-616-P

STAT-BAG™ Heavy-duty Waste Bags

Disposable bags for trash, biohazard waste and linens in the operating room or any area of the facility.

- Infused with specially formulated resin to help prevent tears
- Accommodate irregular shaped items
- Available with drawcord or drawtape for easy, secure closure
- Additional sizes and colors available



ITEM #	PRODUCT NAME	DESCRIPTION
BCD-LG	STAT-BAG™	Clear Bag with Drawtape; Large
BCC-LG	STAT-BAG™	Clear Bag with Drawcord; Large
KBC	STAT-BAG™	Clear Kick Bucket Bag
BRD-LG	STAT-BAG™	Red Bag with Drawtape; Large with Biohazard symbol
BRC-LG	STAT-BAG™	Red Bag with Drawcord; Large with Biohazard symbol
KBR	STAT-BAG™	Red Kick Bucket Bag; with Biohazard symbol
BBD-LG	STAT-BAG™	Blue Bag with Drawtape; Large
BBC-LG	STAT-BAG™	Blue Bag with Drawcord; Large

* Data on file.

Additional Room Turnover Components

A variety of other items are available to meet the specific needs of any department:

- Fast-acting solidifier for safe disposal of blood and body fluids
- Wipes for cleaning surfaces during room turnover
- Urology drain bags designed for a variety of cystoscopy tables (adapters and footswitch covers also available)
- Absorbent floor pads to collect fluid and keep the floor dry to help prevent slips and falls
- Cord covers to secure cords on the floor and help prevent tripping
- Patient transfer sheet designed to reduce friction and help decrease physical strain to the back, neck, shoulders and arms of healthcare workers




ITEM #	PRODUCT NAME	DESCRIPTION
SS-500	STAT-SORB™	Solidifier for 500cc Canister
SS-1200	STAT-SORB™	Solidifier for 1200cc Canister
SS-1500	STAT-SORB™	Solidifier for 1500cc Canister
SS-2000	STAT-SORB™	Solidifier for 2000cc Canister
SS-3000	STAT-SORB™	Solidifier for 3000cc Canister
WDW-L	Wipe	Large, Dry Wipe, White
UDB-10	Urology Drain Bag	Urology Drain Bag for Siemens, Dornier 1000
UDB-20	Urology Drain Bag	Urology Drain Bag for Leibel 2000/OEC 2500
UDB-28	Urology Drain Bag	Urology Drain Bag for GE OEC 2800
UDB-30	Urology Drain Bag	Urology Drain Bag for Skytron 3000
See page 32	Trip-No-More™	Disposable Cord Cover
See page 33	DriFloor™	Fluid Absorption Floor Pad
See page 35	Z-Slider™	Patient Transfer Sheet

Z-Cord Holder

Designed to keep cords and tubes organized on the sterile field.

- Grooved holes allow for easy and secure placement of cords and tubes
- Adhesive backing keeps it firmly in place
- X-ray detectable



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1117	Z-Cord Holder	Cord & Tube Holder (Sterile)	50/Bx 2 Bx/Cs	100

Z-Suture Saddle Bag

A clear plastic bag that allows for easy identification of its contents.

- Can be draped over the back table or Mayo stand
- Perforated so bags can be placed in separate locations



ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1130	Z-Suture Saddle Bag	20" for use on the Mayo Stand or Back Table (Sterile)	50/Bx 2 Bx/Cs	100

Available sterile and non-sterile

SANDEL REFERENCES

1. Jagger J, Berguer R, Phillips EK, Parker G, Gomaa AE. Increase in sharps injuries in surgical settings versus nonsurgical settings after passage of national needlestick legislation. *J Am Coll Surg.* 2010; 210(4): 496-502.
2. OSHA Dealing with Sharps Injuries. <http://www.articlesbase.com/diseases-and-conditions-articles/oshadealing-with-sharps-injuries-5354841.html>. Accessed Jan. 12, 2012.
3. Slip, Trip and Fall Prevention for Healthcare Workers, Dept. of Health and Human Services, Centers for Disease Control and Prevention, <https://www.cdc.gov/niosh/docs/2011-123/pdfs/2011-123.pdf>, Dec. 2010. Accessed 6/9/20.
4. American Nurses Association. 2011. 2011 Health and Safety Survey Report.
5. Safe Patient Handling Programs – Effectiveness and Cost Savings, OSHA. [Osha.gov/dsg/hospitals/documents/3.5_SPH_effectiveness_508.pdf](http://www.osha-slc.gov/dsg/hospitals/documents/3.5_SPH_effectiveness_508.pdf). Accessed 12/2/19.
6. American Nurses Association. <http://www.nursingworld.org/MainMenuCategories/Policy-Advocacy/State/Legislative-Agenda-Reports/State-SafePatientHandling>. Accessed Aug. 25, 2015
7. Van Den Bos J et al. The \$17.1 Billion Proglem: The Annual Cost of Measurable Medical Errors, *Health Affairs.* April 2011 Volume 30 (4): 596-603.
8. The Joint Commission. National Patient Safety Goals Effective January 1, 2015. http://www.jointcommission.org/assets/1/6/2015_NPSG_HAP.pdf. Accessed August 25, 2015
9. Agency for Healthcare Research and Quality. <https://www.ahrq.gov/sites/default/files/publications/files/putoolkit.pdf>. Accessed March 13, 2017.
10. Redistributing weight is a known method for reducing pressure from bony prominences. Overall pressure reduction is known to reduce the risk of pressure ulcers. Test data indicates that overall pressure of ProForm products are lower than that of traditional foam positioning products. Data on file.
11. Data on file, external testing completed by Vascocare Medical Ltd., 2016.
12. "Preventing Healthcare-Associated Infections." CDC. Last accessed 18 Sept 2014; <http://www.cdc.gov/washington/~cdcatWork/pdf/infections.pdf>
13. "Addressing Rising Health Care Costs: A Resource Guide for Consumer Advocates." Robert Wood Johnson Foundation. November 2013; 22.
14. Clemson University, School of Materials Science & Engineering. February 28, 2011. Data report on file.
15. AORN Recommended Practices for Environmental Cleaning in the Perioperative Setting. In: *Perioperative Standards and Recommended Practices*. Denver, CO: AORN, Inc., 2010.





Ansell has a long-standing commitment to advancing education that adds **value beyond product**.

Guided and supported by the collective expertise of leading scientists, physicians and other key opinion leaders worldwide, AnsellCARES research and education initiatives promote safer working and living environments for everyone.



EDUCATION

Industry-leading learning experiences

Accredited and non-accredited continuing education (CE) programs help to advance your professional development as well as encourage best practices and optimal outcomes in healthcare provider and patient safety.

EVIDENCE

Grounded in clinical research

Multimedia resources such as white papers, newsletters, and videos provide current knowledge to help you continue protecting and reducing risks for you and your patients in an ever-changing healthcare environment.



ENGAGEMENT

Connected with healthcare and safety leaders

Our global network of key opinion leaders and leading clinical specialists helps us provide you with timely advice, actionable recommendations, and ongoing education.





Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: info@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0)2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com